

Product Certificate

Item number / Item name

Item No	Item name
654510	FioniaVet Blood transfusion set, 16Gx1½" needle, sterile, 5/pk

Photo



Product type

Dispensing bag 70ml of anticoagulant CPDA-1 solution.

100 ml CPDA-1 contains:

- Citric acid: 0.299g
- Sodium citrate: 2.63g
- Sodium dihydrogen phosphate: 0.222g
- Glucose: 3.19g
- Adenine: 0.0275g
- Water to 100 milliliters.

A 300 milliliter usable blood bag contains 35 milliliters of CPDA-1.

Product purpose

For veterinary use only.

Size

15 * 32 cm

Material type(s)

PVC Bag

Colour

Transparent

Shelf life

3 years

Packing

3-layer sterile packing.

Single packaging and sterile, with 5 single bags packaged in an aluminum bag.

Storage condition

Dry place.

CE marking

No, if prepared for veterinary use only

Recycling and Disposal

Product parts

	Material	Recycling symbol	Percentage of re-used material	Recycling / Disposal
PVC Bag	PVC			

Packaging materials

	Material	Recycling symbol	Percentage of re-used material	Recycling / Disposal
First layer (sterile) Pouch	PET/CPP			Separate materials and send to recycling.
Second layer – dispenser box	PET/AL/PE			The packaging / shipping materials are recyclable, unbleached

Approval

Signed	Function
Martin Eberhardt	e-mail: eberhardt@fioniavet.com

Document History

Version	Date	Change
1	22-may-2025	First release 2025