

PB-1CD456M0Y

1. Intended use

PB-1CD456M0Y is a single blood bag, which contain 63ml CPDA-1 anticoagulant solution and is intended for the collection, storage and transfusion of 450ml human blood. This blood bag system has following safety features: Needle injury protector, Sampling arm, Holder and Pre-donation sampling bag.

2. Packaging

The blood bags are packed in a primary package made of Polypropylene and then in a secondary package of Aluminium foil. These Aluminium foils are then packed in a shipping carton.

Product code	Qty of bags / Aluminum foil	Qty of Aluminium foils / Carton	Qty of bags / Carton
PB-1CD456M0Y	8 units	5 units	40 units

3. Technical Specifications

Specifications

Material	:	PVC (DEHP)
Inner diameter of the tube	:	2.95 ± 0.05 mm
Outer diameter of the tube	:	4.4 ± 0.1 mm
Donor tube length	:	105 ± 5 cm (12 marks)
Phlebotomy needle	:	16-gauge needle

Bag dimensions

Bag Type	Volume (ml)	Inside Width (mm)	Outside Width (mm)	Inside Length (mm)	Outside Length (mm)
Primary bag	450ml	120 ± 3	132 ± 3	180 ± 5	204 ± 5

Drawing

